



neoLab® XY rack red

Product Images



Description

Double-sided rack
with 96 places for
- PCR or 0.5 ml tubes on one side and
- 1.5/2.0 ml tubes on the other side (max. 12 mm Ø).
Sliding tabs on the front and long side in combination with the alphanumeric index facilitate orientation during work.
Graduated recesses give tubes with 1.5 ml more grip.

Material: PP
Temperature resistant to -80°C
Dimensions: 206 x 139 x 27 mm (L x W x H)

No.	2-1959
Manufacturer (Brand)	neoLab
EAN	4058072059057
Color	rot
Material	Polypropylen (PP)
Wide	139 mm
Height	27 mm
suitable for	Mikrogefäße/ Röhrchen
Racks features	stapel-/zusammensteckbar mit besonderer chemikalischer Beständigkeit inkubatorgeeignet für verschiedene Gefäßtypen/Volumina alphanumerisch codiert
Number of floors	1
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Gefrierschrank/ -truhe
Squares (WxL)	8 x 12
Type frame	Gestell für Reaktionsgefäße. PCR-Tubes. Kryoröhrchen
for vessels with DM up to	5.5 mm 11.0 mm
for micro vessel/tube shape	für alle Formen
for type (micro) vessel	Reaktionsgefäße PCR-Reaktionsgefäße

