




## neoLab® XY rack blue

### Product Images

---



### Description

---

Double-sided rack

with 96 places for

- PCR or 0.5 ml tubes on one side and

- 1.5/2.0 ml tubes on the other side (max. 12 mm Ø).

Sliding tabs on the front and long side in combination with the alphanumeric index facilitate orientation during work.

Graduated recesses give tubes with 1.5 ml more grip.

Material: PP

Temperature resistant to -80°C

Dimensions: 206 x 139 x 27 mm (L x W x H)

|                               |                                                                                                                                                    |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| No.                           | 2-1960                                                                                                                                             |
| Manufacturer (Brand)          | neoLab                                                                                                                                             |
| EAN                           | 4058072059064                                                                                                                                      |
| Color                         | blau                                                                                                                                               |
| Material                      | Polypropylen (PP)                                                                                                                                  |
| Wide                          | 139 mm                                                                                                                                             |
| Height                        | 27 mm                                                                                                                                              |
| suitable for                  | Mikrogefäße/ Röhrchen                                                                                                                              |
| Racks features                | stapel-/zusammensteckbar mit besonderer chemikalischer Beständigkeit inkubatorgeeignet für verschiedene Gefäßtypen/Volumina alphanumerisch codiert |
| Number of floors              | 1                                                                                                                                                  |
| for vessels with volume up to | 0.5 ml 1.5 ml 2 ml                                                                                                                                 |
| suitable for cooling in       | Gefrierschrank/ -truhe                                                                                                                             |
| Squares (WxL)                 | 8 x 12                                                                                                                                             |
| Type frame                    | Gestell für Reaktionsgefäße. PCR-Tubes. Kryoröhrchen                                                                                               |
| for vessels with DM up to     | 5.5 mm 11.0 mm                                                                                                                                     |
| for micro vessel/tube shape   | für alle Formen                                                                                                                                    |
| for type (micro) vessel       | Reaktionsgefäße PCR-Reaktionsgefäße                                                                                                                |

