



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® XY rack green

**€32.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Double-sided rack

with 96 places for

- PCR or 0.5 ml tubes on one side and
- 1.5/2.0 ml tubes on the other side (max. 12 mm Ø).

Sliding tabs on the front and long side in combination with the alphanumeric index facilitate orientation during work.

Graduated recesses give tubes with 1.5 ml more grip.

Material: PP

Temperature resistant to -80°C

Dimensions: 206 x 139 x 27 mm (L x W x H)

## Additional Information

No.	2-1961
Manufacturer (Brand)	neoLab
EAN	4058072059071
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	green
Material	Polypropylene (PP)
Wide	139 mm
Height	27 mm
suitable for	Micro Vessels / Tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	8 x 12
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	5.5 mm 11.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Reaction vessels PCR reaction vessels

