



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine reaction tube rack, yellow, PP, 2 x 8 tubes 5 ml

**€18.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

This neoLab reaction tube stand was specially developed for 5 ml reaction tubes. Made of durable polypropylene, this reaction tube stand is temperature resistant from -90°C to +120°C and autoclavable. The feet of the reaction vessel rack are secured with knobs to prevent slipping on the tabletop, but should be removed for autoclaving. With a division of 2 x 8 rows, this neoLab reaction tube rack offers space for 16 x 5 ml reaction tubes or eppis. Conical recesses guarantee a firm hold of the reaction tubes in the rack. Due to the wedge shape, the reaction vessels in the back row can also be easily reached and processed. All openings are numbered.

## Additional Information

---

No.	2-1900
Manufacturer (Brand)	neoLabLine
EAN	4058072058753
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Wide	80 mm
Height	45 mm
Racks features	incubator suitable
Number of floors	1
DM Bore	16.5 mm
Squares (WxL)	2 x 8
Type frame	Rack for reaction tubes. PCR tubes. Cryotubes
for vessels with DM up to	16.0 mm
for vessels with volume up to	5 ml
for micro vessel/tube shape	conical
suitable for cooling in	Refrigerator freezer/ freezer

