




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® water bath insert for 20 cryotubes

€45.00
plus VAT &
Shipping

Product Images



Description

Round water bath insert, floating, for 20 cryotubes of the most common sizes.

The rack is provided with bars which guarantee good stability.

The rack floats in water and is therefore suitable for thawing cryotubes in a water bath.

Dimensions: 6.5 mm thick, 95 mm Ø

Feet: 19 mm long

Material: Polypropylene

Additional Information

No.	1-1740
Manufacturer (Brand)	neoLab
EAN	4058072033385
diameter	95 mm
Color	blue
Material	Polypropylene (PP)
Height	80 mm
DM Hole	12.5 mm
for vessels with volume up to	1.5 ml 2 ml
for vessels with DM up to	12.0 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

