



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Water bath insert, sinking, 3 x 5 R., for 27-30 mm Dchm.

€91.30
plus VAT &
Shipping

Product Images



Description

These water bath inserts sink without adding weight! Molded from a special PP compound, heavier than water. With handle openings on the sides.

Autoclavable up to 121°C.

Additional Information

No.	6-2314
Manufacturer (Brand)	neoLab
EAN	4058072025717
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Temp MAX	121 °C
Temp MIN	0 °C
Wide	11 mm
Height	12.7 mm
Length	21.3 mm
Depth	21.3 mm
DM Hole	31 mm
Capacity	15
Squares (WxL)	3 x 5
for vessels with DM up to	30.0 mm
for vessels with volume up to	30 ml
for type (micro) vessel	Test tubes/tubes

