



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® aluminum block MTP form 2 x 4 hole for 5.0 ml

Special Price

€169.00 was

€210.00

plus VAT &
Shipping

Product Images



Description

Aluminum block for temperature control of reaction tubes and PCR tubes.

- small space requirement
- fits on any ice bath, in common water baths or on heating plates
- optimal for all temperature-critical work
- extremely fast cooling, long cooling effect.
- Use pre-cooled or pre-warmed, even without pre-tempering
- even temperature distribution
- special shape of the holes for optimum heat exchange and avoidance of specimen losses due to falling out
- low-priced special designs on request

Additional Information

No.	7-8004
Manufacturer (Brand)	neoLabLine
EAN	4058072041724
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Aluminum
Wide	80 mm
Form	rectangular
Height	60 mm
Length	115 mm
Depth	115 mm
suitable for	Reaction vessels Micro vessels/ tubes
DM Bore	17.5 mm
Squares (WxL)	2 x 4
for vessels with DM up to	16.0 mm
for vessels with volume up to	5 ml
for type (micro) vessel	Reaction vessels Micro vials

