



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



CellCamper® Alublock 120x0.2 ml PCR


Special Price

€289.00 was

€379.00

plus VAT &
Shipping

Product Images



Description

The CellCamper aluminum blocks are suitable for the CellCamper Maxi and Midi. In general, these aluminum blocks can be used at temperatures from -196 to over 100°C. The temperature distribution over the entire block is very uniform and deviates only by ± 0.1 °C. All aluminum blocks can be easily cleaned with soap suds or alcohols and are stable under UV light.

Capacity: 120x0.2 ml PCR tubes

Additional Information

No.	2-3719
Manufacturer (Brand)	CellCamper
EAN	4058072090074
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	dark gray
Material	Aluminum
Wide	94 mm
Form	rectangular
Height	17 mm
Length	135 mm
Depth	135 mm
suitable for	Reaction vessels
DM Bore	6 mm
for vessels with volume up to	0.2 ml
Squares (WxL)	10 x 12
for vessels with DM up to	5.5 mm
for type (micro) vessel	PCR reaction vessels

