



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



CellCamper® Alublock 35x1.5 ml reaction tubes


Special Price

€219.00 was

€290.00

**plus VAT &
Shipping**

Product Images



Description

The CellCamper aluminum blocks are suitable for the CellCamper Maxi and Midi. In general, these aluminum blocks can be used at temperatures from -196 to over 100°C. The temperature distribution over the entire block is very uniform and deviates only by ± 0.1 °C. All aluminum blocks can be easily cleaned with soap suds or alcohols and are stable under UV light.

Capacity: 35x1.5 ml reaction tubes

Additional Information

No.	2-3720
Manufacturer (Brand)	CellCamper
EAN	4058072090081
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	dark gray
Material	Aluminum
Wide	94 mm
Form	rectangular
Height	38 mm
Length	135 mm
Depth	135 mm
suitable for	Reaction vessels
DM Bore	12 mm
for vessels with volume up to	1.5 ml
Squares (WxL)	5 x 7
for vessels with DM up to	11.5 mm
for type (micro) vessel	Reaction vessels

