



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## CellCamper® Alublock 35x2.0 ml reaction tubes

Special Price

**€219.00** was

**€290.00**

plus VAT &  
Shipping

### Product Images



## Description

---

The CellCamper aluminum blocks are suitable for the CellCamper Maxi and Midi. In general, these aluminum blocks can be used at temperatures from -196 to over 100°C. The temperature distribution over the entire block is very uniform and deviates only by  $\pm 0.1$  °C. All aluminum blocks can be easily cleaned with soap suds or alcohols and are stable under UV light.

**Capacity: 35x2.0 ml reaction tubes**

## Additional Information

---

No.	2-3721
Manufacturer (Brand)	CellCamper
EAN	4058072090098
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	dark gray
Material	Aluminum
Wide	94 mm
Form	rectangular
Height	38 mm
Length	135 mm
Depth	135 mm
suitable for	Reaction vessels
DM Bore	12.5 mm
Capacity	35
Squares (WxL)	5 x 7
for vessels with DM up to	12.5 mm
for vessels with volume up to	2 ml
for type (micro) vessel	Reaction vessels

