



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



CellCamper® Alublock 35x2.0 ml Cryovials

Special Price

€266.00 was

€349.00

**plus VAT &
Shipping**

Product Images



Description

The CellCamper aluminum blocks are suitable for the CellCamper Maxi and Midi. In general, these aluminum blocks can be used at temperatures from -196 to over 100°C. The temperature distribution over the entire block is very uniform and deviates only by ± 0.1 °C. All aluminum blocks can be easily cleaned with soap suds or alcohols and are stable under UV light.

Capacity: 35x2.0 ml cryovials

Additional Information

| | |
|-------------------------------|----------------------------------|
| No. | 2-3722 |
| Manufacturer (Brand) | CellCamper |
| EAN | 4058072090104 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | dark gray |
| Material | Aluminum |
| Wide | 94 mm |
| Form | rectangular |
| Height | 38 mm |
| Length | 135 mm |
| Depth | 135 mm |
| suitable for | Reaction vessels Cryogenic tubes |
| DM Bore | 12.5 mm |
| for vessels with volume up to | 2 ml |
| Squares (WxL) | 5 x 7 |
| for vessels with DM up to | 12.5 mm |
| for type (micro) vessel | Cryotube |

