



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® reaction block, ¼ circle segment, 4 hole

€59.20
plus VAT &
Shipping

Product Images



Description

Anodized aluminum blocks with different hole diameters for heating reaction vessels, ¼ cake piece each. The blocks are suitable for placing on hot plates and magnetic stirrers. Corresponding glass reaction vessels available on request.

Additional Information

No.	D-8304
Manufacturer (Brand)	sunlab
EAN	4058072024314
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Aluminum anodized
Form	1/4 piece of cake
suitable for	Test tubes Centrifuge tubes
DM Bore	28 mm
for vessels with volume up to	40 ml
for vessels with DM up to	28.0 mm
for type (micro) vessel	Test tubes/tubes

