



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf SmartBlock PCR 384, Thermoblock for PCR Plates 384, incl. Lid

**€805.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Eppendorf SmartBlock™ PCR 384, Thermoblock for PCR Plates 384, incl. Lid. Lid - For 384-well PCR plates - Optimized block design ensures optimal temperature transfer within your sample for excellent results - Integrated sensor technology for Eppendorf SmartBlock-specific calibration for highest temperature accuracy and homogeneity - Each SmartBlock is individually tested for temperature accuracy during production and delivered with a serial number-specific certificate for documentation purposes - Insulation of the thermoblocks for ergonomic working - without burning fingers (PhysioCare)specific certificate - Insulation of the thermoblocks for ergonomic working - without burning the fingers (PhysioCare Concept) - Automatic block recognition for safe handling - Eppendorf QuickRelease™ technology enables fast and easy exchange of the blocks - without any tools Only usable for ThermoMixer C and ThermoStat C

## Additional Information

---

No.	VB-1780
Manufacturer (Brand)	Eppendorf
EAN	4058072170752
Transport temperature	Room temperature
Color	gray
Material	Aluminum
Racks features	with closure/lid
Number of floors	1

