



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf SmartBlock™ 0.5 ml, thermoblock for 24x 0.5 ml

€648.00
plus VAT &
Shipping

Product Images



Description

Eppendorf SmartBlocks

for all common tube formats in the volume range from 5 µl to 50 ml.

Also for PCR plates.

The block can be changed quickly by means of a simple mechanism.

A plastic sheath provides thermal insulation and protection against accidental contact.

Order No. E-5055 and E-5056: with lid

More SmartBlocks in the OnlineShop and on request.

Additional Information

No.	VB-1795
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-5050
EAN	4058072418557
Transport temperature	Room temperature
Color	gray
Material	Aluminum
Form	rectangular
suitable for	Reaction vessels
Number of floors	1
DM Bore	6.2 mm
Material block	Aluminum
for vessels with volume up to	0.5 ml
Squares (WxL)	4 x 6
for vessels with DM up to	8.0 mm
for type (micro) vessel	Reaction vessels

