



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf SmartBlock™ 2.0 ml, thermoblock for 24x 2.0 ml

**€648.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Eppendorf SmartBlocks

for all common tube formats in the volume range from 5 µl to 50 ml.

Also for PCR plates.

The block can be changed quickly by means of a simple mechanism.

A plastic sheath provides thermal insulation and protection against accidental contact.

Order No. E-5055 and E-5056: with lid

More SmartBlocks in the OnlineShop and on request.

## Additional Information

---

No.	VB-1796
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-5052
EAN	4058072168339
Transport temperature	Room temperature
Color	gray
Material	Aluminum
Form	rectangular
suitable for	Reaction vessels
Number of floors	1
DM Bore	11.2 mm
Capacity	24
Material block	Aluminum
Squares (WxL)	4 x 6
for vessels with DM up to	11.0 mm
for vessels with volume up to	2 ml
for type (micro) vessel	Reaction vessels

