



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine reaction exchange blocks WB for 15 x 0.5 and 20 x 1.5 ml reaction tubes

€295.00
plus VAT &
Shipping

Product Images



Description

Reaction exchange block WB 0.5/1.5 ml

For 15 x 0.5 and 20 x 1.5 ml reaction tubes

Additional Information

| | |
|-------------------------------|------------------|
| No. | 7-0926 |
| Manufacturer (Brand) | neoLabLine |
| EAN | 4058072158644 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Speed MAX | 1500 rpm |
| Color | gray |
| Wide | 280 mm |
| Height | 170 mm |
| Length | 220 mm |
| Depth | 220 mm |
| suitable for | Reaction vessels |
| Racks features | with closure/lid |
| Material block | Aluminum |
| Squares (WxL) | 5 x 7 |
| for vessels with DM up to | 8.0 mm 11.0 mm |
| for vessels with volume up to | 0.5 ml 1.5 ml |

