



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rack for centrifuge tubes 15 ml, 15 places

€45.10
plus VAT &
Shipping

Product Images



Description

Stable wire rack with epoxy coating for 15 centrifuge tubes with 15 ml capacity.

The additional grids at the bottom guarantee a safe stand especially for conical tubes, but also for all other types and sizes of tubes up to 19 mm diameter.

Additional Information

No.	2-6618
Manufacturer (Brand)	neoLab
EAN	4058072008147
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Steelepoxyc coated
Wide	76 mm
Height	76 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
Tray size	19 x 19 mm
Capacity	15
Squares (WxL)	3 x 5
Type frame	Centrifuge tube rack
for vessels with DM up to	19.0 mm
for vessels with volume up to	15 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

