



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Moonlab® Combi rack for centrifuge tubes,  
foldable, 15 ml: 21 places (7x3) / 50 ml: 12 places  
(6x2), blue, PP, 1 pc**

**€42.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

This foldable rack for centrifuge tubes is practical and can be stored in a space-saving way. All tubes with a diameter of 16 or 30 mm can be organized here. By simply unfolding the rack stands stable. With handles on the sides and alphanumeric index, for easier identification of the tubes.

The rack is made of sturdy polypropylene and is autoclavable. Available in blue color.

## Additional Information

---

No.	4-0064
Manufacturer (Brand)	Moonlab
EAN	4058072085223
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polypropylene (PP)
Wide	137 mm
Height	72 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	16 mm
Capacity	2 12
Type frame	Centrifuge tube rack
for vessels with DM up to	17.0 mm 30.0 mm
for vessels with volume up to	15 ml 50 ml
suitable for cooling in	Refrigerator

