



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




**Moonlab® Combi rack for centrifuge tubes,
foldable, 15 ml: 21 places (7x3) / 50 ml: 12 places
(6x2), pink, PP, 1 pc**

**Price on
request**

**plus VAT &
Shipping**

Product Images



Description

This foldable rack for centrifuge tubes is practical and can be stored in a space-saving way. All tubes with a diameter of 16 or 30 mm can be organized here. By simply unfolding the rack stands stable. With handles on the sides and alphanumeric index, for easier identification of the tubes.

The rack is made of sturdy polypropylene and is autoclavable. Available in the color pink.

Additional Information

No.	4-0065
Manufacturer (Brand)	Moonlab
EAN	4058072085230
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	pink
Material	Polypropylene (PP)
Wide	137 mm
Height	72 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	16 mm
for vessels with volume up to	15 ml 50 ml
suitable for cooling in	Refrigerator
Type frame	Centrifuge tube rack
for vessels with DM up to	17.0 mm 30.0 mm

