



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




Moonlab® Combi rack for centrifuge tubes, ABS, 30 x 15 ml, 20 x 50 ml

€8.50 As low as

€7.20

**plus VAT &
Shipping**

Product Images



Description

This **rack** is suitable for storing conical 15 or 50 ml centrifuge tubes. Constructed of ABS plastic for stable use in a water bath. The holes are alphanumerically marked.

Additional Information

No.	4-0053
Manufacturer (Brand)	Moonlab
EAN	4058072085193
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	orange
Material	Acrylonitrile butadiene styrene (ABS)
sterile	No
Wide	171 mm
Height	63 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	30 mm
Squares (WxL)	4 x 5 5 x 6
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	17.0 mm 30.0 mm
for vessels with volume up to	15 ml 50 ml
for micro vessel/tube shape	conical round
suitable for cooling in	Refrigerator freezer/ freezer

