



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® Combi rack for centrifuge tubes, PP, 18 x 15 ml, 10 x 50 ml

€6.00
plus VAT &
Shipping

Product Images



Description

The perfect **storage system** for conical **centrifuge tubes** with a capacity of 18x 15 ml and 10 x 50 ml. Made of durable PP material, this combination rack provides a sturdy and reliable solution for organizing your samples. The **alphanumeric markings** on the holes make it easy to identify and manage your samples. Order the Moonlab® centrifuge tube rack now and optimize your lab work.

Additional Information

No.	4-0092
Manufacturer (Brand)	Moonlab
EAN	4058072118600
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
Wide	102 mm
Height	63 mm
Racks features	for different vessel types/volumes alphanumerically coded
Number of floors	1
Squares (WxL)	3 x 9 2 x 5
Type frame	Combi rack for different vessels and volumes
for vessels with DM up to	17.0 mm 30.0 mm
for vessels with volume up to	15 ml 50 ml
for micro vessel/tube shape	for all shapes
suitable for cooling in	Refrigerator freezer/ freezer

