



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Pop-Up™ Rack for 15 and 50 ml, pink, 2 pcs./pack

€49.00
plus VAT &
Shipping

Product Images



Description

Stored flat, unfolded quickly!

When the rack is not in use, it folds up very small. When it is needed, it is unfolded with a flick of the wrist and holds cell culture tubes of 15 or 50 ml.

The Pop-Up® Rack holds 21 tubes of 15 ml and 12 tubes of 50 ml, each space is marked separately.

The rack is stackable in folded and set-up state.

Material PP, not recommended for autoclaving

Dimensions: 255 x 137 mm (L x D),
height folded 21 mm, erected 72 mm

Additional Information

No.	2-6625
Manufacturer (Brand)	neoLab
EAN	4058072008215
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	pink
Material	Polypropylene (PP)
Wide	137 mm
Height	72 mm
suitable for	Test tubes/tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes alphanumerically coded
Number of floors	1
DM Bore	30 mm
for vessels with volume up to	15 ml 50 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	variable
Type frame	Centrifuge tube rack
for vessels with DM up to	19.0 mm 30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Centrifuge tube

