



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Acrylic stand for conical 50 ml tubes, 2 x 6 hole

€150.10
plus VAT &
Shipping

Product Images



Description

Stand made of acrylic glass, heavy duty, glued and screwed.
Suitable for conical 50 ml Falcon centrifuge tubes 2070 or "Blue Max".

Hole diameter top and center: 30 mm

Hole diameter bottom: 20 mm

Height: 85 mm

Additional Information

No.	2-2281
Manufacturer (Brand)	neoLab
EAN	4058072060305
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	85 mm
Height	85 mm
suitable for	Test tubes/tubes
Racks features	stackable/connectable for different vessel types/volumes
Number of floors	1
DM Bore	30 mm
for vessels with volume up to	50 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	2 x 6
Type frame	Centrifuge tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	conical
for type (micro) vessel	Centrifuge tube

