



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Stand for centrifuge tubes 50 ml, 2 x 4 hole

€53.40
plus VAT &
Shipping

Product Images



Description

Wire rack for conical 50 ml "Blue Max" centrifuge tubes, specially epoxy coated.

Autoclavable.

Tray size: 30 x 30 mm

8 openings (2 x 4):

147 x 80 x 89 mm high

Additional Information

No.	1-3110
Manufacturer (Brand)	neoLab
EAN	4058072035044
autoclavable at 121°C	Yes
Color	blue
Material	Steel epoxy coated
Wide	80 mm
Height	89 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
Tray size	30 x 30 mm
Squares (WxL)	2 x 4
Type frame	Centrifuge tube rack
for vessels with DM up to	30.0 mm
for vessels with volume up to	50 ml
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

