



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® centrifuge bottle rack for 6 bottles

**€49.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For storing conical centrifuge bottles with a capacity of 250 ml (e.g. Falcon, Corning, Nunc).

Made of steel wire covered with synthetic resin.

The 60 mm square grids provide space for 6 vessels.

The lower level has rails for a better hold of the bottle bottoms.

Dimensions: 208 x 135 x 102 mm

## Additional Information

No.	6-2130
Manufacturer (Brand)	neoLab
EAN	4058072025335
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Steelepoxyc coated
Wide	135 mm
Height	102 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
for vessels with volume up to	250 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	60 x 60 mm
Squares (WxL)	2 x 3
Type frame	Centrifuge tube rack
for vessels with DM up to	60.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Centrifuge bottles

