



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Centrifuge tube holder f. 50 ml, 5 pcs./set

€45.30
plus VAT &
Shipping

Product Images



Description

5 interlocking holders for 50 ml centrifuge tubes in a set.

Each holder can be used individually.

Round and conical tubes are held securely in the shaped base recess.

Autoclavable.

The set consists of one holder each in blue, green, yellow, pink and orange.

Dimensions: 53 x 70 x 67 mm

Additional Information

No.	2-6612
Manufacturer (Brand)	neoLab
EAN	4058072008116
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue green yellow pink orange
Material	Polypropylene (PP)
Wide	53 mm
Height	70 mm
suitable for	Test tubes/tubes
Racks features	stackable/assemblable with special chemical resistance suitable for incubator for different vessel types/volumes
Number of floors	1
DM Bore	30 mm
for vessels with volume up to	50 ml
suitable for cooling in	Freezer/ freezer
Squares (WxL)	5 x 1
Type frame	Centrifuge tube rack
for vessels with DM up to	30.0 mm
for micro vessel/tube shape	for all shapes
for type (micro) vessel	Centrifuge tube

