



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Stand for 1 centrifuge tube 50 ml

€21.70
plus VAT &
Shipping

Product Images



Description

Stand for a single 50 ml centrifuge tube.

Material:

Wire, coated with epoxy resin.

Easy handling and safe stand due to nest for tube bottom.

Total height: 67 mm

Bottom Ø: 67 mm

Additional Information

No.	1-3108
Manufacturer (Brand)	neoLab
EAN	4058072035037
autoclavable at 121°C	Yes
Color	orange
Material	Steel epoxy coated
Wide	67 mm
Height	67 mm
suitable for	Test tubes/tubes
Racks features	with special chemical resistance suitable for incubators for different vessel types/volumes
Number of floors	1
DM Bore	30 mm
Squares (WxL)	1 x 1
Type frame	Centrifuge tube rack
for vessels with DM up to	30.0 mm
for vessels with volume up to	50 ml
for micro vessel/tube shape	conical
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

