



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Benchtop Stand for Autoclave Bags, 125 x 75 x 205 mm

€45.60
plus VAT &
Shipping

Product Images



Description

Table stand made of steel wire

epoxy resin coated,
impact resistant and highly chemically resistant.
Suitable for bags 22 x 28 cm.
Dimensions: 125 x 75 x 205 mm (W x D x H)

Accessories for Biohazard autoclave bags (Order No. 6-1042)

Additional Information

No.	6-1043
Manufacturer (Brand)	neoLab
EAN	4058072024697
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	orange
Material	Steel epoxy coated
Volume	1500 ml
Wide	75 mm
Height	205 mm
Length	125 mm
Depth	125 mm
suitable for	variable applications

