



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® four feet for glass ceramic laboratory protection plates made of chrome nickel steel 155 x 155 mm

€58.19 As low as
€46.55
plus VAT &
Shipping

Product Images



Description

Square base, for plate 155 x 155 mm In the laboratory, solutions are often heated. This step can take place using various techniques, e.g. in a water or sand bath, with the aid of a heating hood or classically using a Bunsen burner. If the Bunsen burner is used as the heating source, it is recommended to use a stand or square base and an appropriate glass ceramic plate. The vessel to be heated is placed on the glass ceramic plate, the Bunsen burner underneath. Dirt on the glass ceramic plate can be removed with the help of the cleaning scraper. Special features and advantages: Holder plate: steel 1.4301, stand feet: steel 1.4305 Made of temperature-resistant chrome-nickel steel Very good stability due to four feet and stable design Application: Plate holder for glass ceramic plates, Heating, For glass ceramic laboratory protection plate.

Additional Information

No.	SD-2599
Manufacturer (Brand)	DURAN
VGKL number	322907757
EAN	4058072184865
Manufacturer EAN	4032051061583
Gross weight	0.37
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.32000 kg
Country of origin	Germany
sterile	No
Region of origin	Rhineland-Palatinate
Packaging volume	0.01530900 cc
Customs tariff number	73269085

