



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Stand for flasks with 50 ml, 3 pcs./pack

**€39.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Stand made of plastic for small round, pointed or pear-shaped flasks.  
The base is slightly conical in shape and has a hole at the bottom for a good fit.  
High stability due to a wide base.

## Additional Information

---

No.	6-5002
Manufacturer (Brand)	neoLab
EAN	4058072026431
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	70 mm
Color	red
Material	Polypropylene (PP)
Height	75 mm
suitable for	Piston
for vessels with volume up to	50 ml
suitable for cooling in	Freezer/ freezer
Type frame	Piston rack
for piston type	Round-bottom Spherical Bulb-shaped pistons

