



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Stand for flasks with 100 ml, 3
pcs./pack**

**€38.80
plus VAT &
Shipping**

Product Images



Description

Stand made of plastic for small round, pointed or pear-shaped flasks.
The base is slightly conical in shape and has a hole at the bottom for a good fit.
High stability due to a wide base.

Additional Information

No.	6-5003
Manufacturer (Brand)	neoLab
EAN	4058072026448
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	80 mm
Color	red
Material	Polypropylene (PP)
Height	90 mm
suitable for	Piston
Type frame	Piston rack
for vessels with volume up to	100 ml
for piston type	Round-bottom Spherical Bulb-shaped pistons
suitable for cooling in	Freezer/ freezer

