




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® electrophoresis plate stand

€133.00
plus VAT &
Shipping

Product Images



Description

Electrophoresis plate stand made of steel with epoxy resin coating for storage and rapid drying of 10 electrophoresis plates.

The support is covered with a rubber sleeve to prevent chipping at the edge of the plate. Rubber buffers on the underside of the stand prevent unintentional slipping on the stand surface.

Dimensions: 20 x 16 x 18 cm high

Suitable for plates with a size of 16.5 cm upwards.

Additional Information

No.	7-9730
Manufacturer (Brand)	neoLab
EAN	4058072050191
Transport temperature	Room temperature
Color	light blue
Material	Steelepoxy coated
Wide	160 mm
Height	180 mm
Length	200 mm
Depth	200 mm
Type frame	Electrophoresis plate rack

