




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Maxi gel drying frame complete in set, 337x250x4 mm

Special Price

€65.00 was

€97.50

**plus VAT &
Shipping**

Product Images



Description

Easy drying of gels by means of frame and between 2 compostable biodegradable cellophane films. The films are permeable and uncoated. The simple and inexpensive system allows drying gels without heat or vacuum. Wrap the gel in foil and fix it with 4 clamps between the plexiglass frames.

Scope of delivery: 1 closed base plate 4 mm thick, 1 frame Maxi DIN A4, 4 mm thick; 4 clamps, 1 pack á 50 sheets of cellophane foil 240 mm wide, 327 mm long

MAXI gel drying frame: Outer dimension: 337 x 250 x 4 mm Inner dimension: 297 x 210 x 4 mm

Additional Information

No.	7-9740
Manufacturer (Brand)	neolab
EAN	4058072050207
Transport temperature	Room temperature
Material	Soda-lime glass (normal glass) 1
Wide	4 mm
Width/inside	4 mm
Height	250 mm
Height inside	210.0 mm
Length	337 mm
Length inside	297 mm
Thickness	4 mm
Spare part	2-4028 7-9743

