



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Midi gel drying frame, complete in set, 250x188x4 mm

€16.24
plus VAT &
Shipping

Product Images



Description

Easy drying of gels by means of frame and between 2 compostable biodegradable cellophane films. The films are permeable and uncoated. The simple and inexpensive system allows drying gels without heat or vacuum. Wrap the gel in foil and fix it with 4 clamps between the plexiglass frames.

Scope of delivery: 1 closed base plate 4 mm thick, 1 frame Midi DIN A5, 4 mm thick, 4 clamps, 1 pack á 50 sheets cellophane foil 178 mm wide, 240 mm long

MIDI gel drying frame: Outer dimension: 250 x 188 x 4 mm Inner dimension: 210 x 150 x 4 mm

Additional Information

No.	7-9741
Manufacturer (Brand)	neolab
EAN	4058072050214
Transport temperature	Room temperature
Material	Soda-lime glass (normal glass) 1
Wide	4 mm
Width/inside	4 mm
Height	188 mm
Height inside	150.0 mm
Length	250 mm
Length inside	210 mm
Thickness	4 mm
Spare part	2-4028 7-9744

