



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Mini gel drying frame, complete in set, 188x145x4 mm

€13.52
plus VAT &
Shipping

Product Images



Description

Easy drying of gels by means of frame and between 2 compostable biodegradable cellophane films. The films are permeable and uncoated. The simple and inexpensive system allows drying gels without heat or vacuum. Wrap the gel in foil and fix it with 4 clamps between the plexiglass frames.

Scope of delivery: 1 closed base plate 4 mm thick ,1 frame Mini DIN A6, 4 mm thick, 4 clamps, 1 pack á 50 sheets of cellophane foil 135 mm wide, 178 mm long.

MINI gel drying frame: Outer dimension: 188 x 145 x 4 mm Inner dimension: 150 x 105 x 4 mm

Additional Information

No.	7-9742
Manufacturer (Brand)	neolab
EAN	4058072050221
Transport temperature	Room temperature
Material	Soda-lime glass (normal glass) 1
Wide	4 mm
Width/inside	4 mm
Height	145 mm
Height inside	104.0 mm
Length	188 mm
Length inside	150 mm
Thickness	4 mm
Spare part	2-4028 7-9745

