



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand test tube rack, PP, for 55 tubes up to $\varnothing$ 18 mm, red

**€67.25**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Sturdy construction, stackable. Very versatile and inexpensive rack for reagent tubes in all common diameters. Filled PP, density 1.2 g/cm<sup>3</sup>, therefore no floating in the water bath. Alphanumerically marked positions for clear identification of samples. Increased spacing between the individual positions. Service temperature -20 °C to 90 °C. Autoclavable (121 °C). The racks are firmly and inseparably mounted in just a few steps.

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | BN-1259          |
| Manufacturer (Brand)  | Brand            |
| VGKL number           | 294943412        |
| EAN                   | 4058072206376    |
| Gross weight          | 1.1225           |
| HSNumber              | 39269097         |
| Storage temperature   | Room temperature |
| Transport temperature | Room temperature |
| Net weight            | 1.09850 kg       |
| Country of origin     | Germany          |
| Packing width         | 0.28000 m        |
| Packing height        | 0.06000 m        |
| Packing depth         | 0.40500 m        |
| Packaging volume      | 6804.000000 ccm  |
| Customs tariff number | 39269097         |

