



KERN stereomicroscope stand OZB,

Product Images



Description

- With our universal stands and basic stands, as well as the microscope heads and external illumination units, you can configure your microscope yourself according to your ideas and tailored to your application
- Thanks to the flexibly applicable and adaptable universal stands, optimal work is possible in all areas with the most varied specimens
- The large universal stands are available as a stand variant and optionally with a table clamp for the edge of the table or for central attachment to your desk. Depending on the model, you can choose between a telescopic arm, an articulated arm or a ball-bearing double-arm universal stand
- A spring-loaded universal stand incl. table clamp simplifies your daily work with your stereomicroscope

Additional Information

| | |
|----------------------|---------------|
| No. | KP-5211 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072147501 |
| Manufacturer | KERN |

