



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN holder for stereomicroscope stand with adjustable handwheel rotation strength

€180.00
plus VAT &
Shipping

Product Images



Description

- There are two microscope head mounts to choose from for these flexible, modular systems. These holders are suitable for all stereomicroscope and universal stands (except spring-loaded arm) to allow precise focusing
- The first option is a holder with adjustable handwheel as well as the adjustment of the rotation strength for your configuration
- For professional applications you can choose a holder with coarse and fine adjustment for optimal focusing
- Diameter of the connection for the microscope head: 76 mm
- Diameter of the connection for the stand column: 25 mm

Additional Information

No.	KP-5215
Manufacturer (Brand)	KERN
EAN	4058072147525
Transport temperature	Room temperature
Manufacturer	CORE

