



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Single stand for 1 pipette Transferpette® electronic 0.25-5 ml or Transferpette® S 2 ml and 0.5-5 ml

€51.05
plus VAT &
Shipping

Product Images



Description

With the stands and shelf holders, you can position your single- or multichannel pipette from BRAND safely and in a space-saving manner in the laboratory. You will find the right solution for every model and every requirement at BRAND:

-

Shelf and wall holders are available for the mechanical single- and multichannel pipettes Transferpette® S.-

The single stands for the Transferpette® S are mostly used for multichannel pipettes.-

For the Transferpette® electronic pipettes, individual stands are available for both single-channel and multichannel pipettes.-

The table stand is a particularly space-saving way to store up to six Transferpette® S or Transferpette® S -8/-12 pipettes.-

Up to three devices (up to 1000 µl) can be charged simultaneously in the charging stand for the Transferpette® electronic. Individual power supply units can then be dispensed with.

Additional Information

No.	BN-1798
Manufacturer (Brand)	Brand
VGKL number	104975386
EAN	4058072210595
Gross weight	0.2055
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.16950 kg
Country of origin	Germany
Packing width	0.11700 m
Packing height	0.08100 m
Packing depth	0.13900 m
Packaging volume	1317.30000000 ccm
Customs tariff number	39269097

