



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand Single rack for 1 Transferpette® -8/ -12 each, all versions

**€25.85**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With the stands and shelf holders, you can position your single- or multichannel pipette from BRAND safely and in a space-saving manner in the laboratory. You will find the right solution for every model and every requirement at BRAND:

-

Shelf and wall holders are available for the mechanical single- and multichannel pipettes Transferpette® S.-

The single stands for the Transferpette® S are mostly used for multichannel pipettes.-

For the Transferpette® electronic pipettes, individual stands are available for both single-channel and multichannel pipettes.-

The table stand is a particularly space-saving way to store up to six Transferpette® S or Transferpette® S -8/-12 pipettes.-

Up to three devices (up to 1000 µl) can be charged simultaneously in the charging stand for the Transferpette® electronic. Individual power supply units can then be dispensed with.

## Additional Information

---

No.	BN-1638
Manufacturer (Brand)	Brand
EAN	4058072209544
Gross weight	0.2324
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.18220 kg
Country of origin	Germany
Packing width	0.13500 m
Packing height	0.08200 m
Packing depth	0.13500 m
Packaging volume	1494.450000 ccm
Customs tariff number	39269097

