



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Charging Stand 2 for an electronic multipette (E3/E3x or stream/Xstream)

€275.00
plus VAT &
Shipping

Product Images



Description

Protect your precious liquid handling instrument from possible damage, spills or chemicals spilled on the work surface. Single instrument stands have a robust design and safely hold an electronic device. With the Charging Stand 2, your Eppendorf Multipette E3/E3x or Multipette stream/Xstream are always fully charged and ready for use. Operate with the plug-in power supply included with the pipette.

Additional Information

No.	VB-0737
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-1973
EAN	4058072409166
Transport temperature	Room temperature
Color	light gray
Height	350 mm
for model pipette	Multipette E3 Multipette E3x
for pipettes / dispensers from	Eppendorf
Type frame	Microliter pipette rack

