



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand 3 stand incl. power supply

**€201.40
plus VAT &
Shipping**

Product Images



Description

3 pc. stand for ^{Transferpette®} electronic (up to 1000 µl) with power supply.

Accessories for ^{Transferpette®} electronic (Order No. 3-0851 ff.)

Additional Information

No.	BN-1799
Manufacturer (Brand)	Brand
old neoLab article no.	3-0857
EAN	4058072210601
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	gray white
suitable for	Pipettes
for model pipette	Transferpette S type variable. 20µl Transferpette S type variable. 100µl Transferpette S type variable. 200µl Transferpette S type variable. 1000µl Transferpette S type variable. 5000µl Transferpette S type variable. 10000µl
for pipettes / dispensers from	Fire
Type frame	Pipette carousel/ loading stand

