



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pipette rack (PP) rectangular, for 50 pipettes

€58.90
plus VAT &
Shipping

Product Images



Description

Space-saving pipette storage due to the vertical arrangement on 5 x 10 spaces.

Openings of the upper and middle decks each have 16 mm Ø.

No contact of the pipette tips with the table top.

Side plates with handles.

Dimensions: 213 x 114 x 222 mm high

Material: Polypropylene

Autoclavable up to 121°C

Additional Information

No.	6-2275
Manufacturer (Brand)	neoLab
EAN	4058072025625
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Wide	213 mm
Height	222 mm
Length	114 mm
Depth	114 mm
suitable for	Pipettes
DM Bore	16 mm
Type frame	Pipette rack

