



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Rack for pipettes 0.1 to 25 ml, acrylic transparent

**€193.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Sturdy pipette rack with 4 shelves for pipettes from 0.1 to 25 ml.  
The acrylic rack is transparent and can be viewed from each side.

Dimensions: 88 x 291 x 401 mm (W x D x H)

## Additional Information

---

No.	1-1000
Manufacturer (Brand)	neoLab
EAN	4058072069513
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Wide	88 mm
Height	401 mm
Length	291 mm
Depth	291 mm
suitable for	Pipettes
Type frame	Pipette rack

