



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rack for pipettes 0.1 to 25 ml, ABS plastic blue

€127.10
plus VAT &
Shipping



Description

Sturdy pipette rack with 4 shelves for pipettes from 0.1 to 25 ml.

The blue rack made of ABS plastic has a cut-out in each shelf for a better view and also a removable holder on each side for storing pipetting balls and pipetting aids within easy reach.

Dimensions: 88 x 291 x 401 mm (W x D x H)

Additional Information

No.	1-1001
Manufacturer (Brand)	neoLab
EAN	4058072069520
autoclavable at 121°C	No
Color	blue
Material	Acrylonitrile butadiene styrene (ABS)
Wide	88 mm
Height	401 mm
Length	291 mm
Depth	291 mm
suitable for	Pipettes
Capacity	variable (depending on vessel volume)
Type frame	Pipette rack

