



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex In-line universal pipette holder 337, for 7 microliter pipettes, pastel pink

€100.78
plus VAT &
Shipping

Product Images



Description

Pipette holder In-line universal 337

Attractive space-saving pipette stand holds up to seven single-channel pipettes of most brands. Made of sturdy polyamide, easy to clean. Stable stand with anti-slip plates.

Available in five pastel colors: ice blue, light gray, mint green, pastel pink, vanilla yellow.

Additional Information

No.	SR-0405
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7848
EAN	4058072369187
Transport temperature	Room temperature
Color	pink
Material	Polyamide (PA)
Type frame	Pipette rack
for pipettes / dispensers from	Socorex Eppendorf Gilson Brand

