



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex In-line universal pipette holder 337, for 7 microliter pipettes, vanilla yellow

**€100.78**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Pipette holder In-line universal 337

Attractive space-saving pipette stand holds up to seven single-channel pipettes of most brands. Made of sturdy polyamide, easy to clean. Stable stand with anti-slip plates.

Available in five pastel colors: ice blue, light gray, mint green, pastel pink, vanilla yellow.

## Additional Information

---

No.	SR-0406
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7849
EAN	4058072369194
Transport temperature	Room temperature
Color	yellow
Material	Polyamide (PA)
for pipettes / dispensers from	Socorex Eppendorf Gilson Brand
Type frame	Pipette rack

