



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Twister™ universal 336, rotating pipette stand for 6 microliter pipettes, emerald green

**€123.97**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Socorex **pipette holder Twister™universal 336**

Pipette stand with unique design and advantageous features.

For six pipettes of most brands, it offers an ideal storage concept for any laboratory. Smooth 360° axis rotation allows quick access to any instrument. Easy disassembly and maintenance.

## Additional Information

---

No.	SR-0396
Manufacturer (Brand)	Socorex
old neolab article no.	7-7839
EAN	4058072369125
Transport temperature	Room temperature
Color	green
Material	Polyamide (PA)
Type frame	Pipette rack
for pipettes / dispensers from	Socorex Eppendorf Gilson Brand

