



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Twister™ universal 336, rotating pipette stand for 6 microliter pipettes, citrin yellow

€123.97
plus VAT &
Shipping

Product Images



Description

Socorex **pipette holder Twister™universal 336**

Pipette stand with unique design and advantageous features.

For six pipettes of most brands, it offers an ideal storage concept for any laboratory. Smooth 360° axis rotation allows quick access to any instrument. Easy disassembly and maintenance.

Additional Information

No.	SR-0401
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7844
EAN	4058072163167
Transport temperature	Room temperature
Color	yellow
Material	Polyamide (PA)
Type frame	Pipette rack
for pipettes / dispensers from	Socorex Eppendorf Gilson Brand

