



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex shelf pipette holder 332 for 2 Socorex microliter pipettes, 4/pack

€67.81
plus VAT &
Shipping

Product Images



Description

Shelf pipette holder 332

Extremely attractive design holds two micro, macro or multichannel pipettes of all Acura® models.

Self-adhesive strip ensures lasting adhesion to any surface.

Made of polycarbonate; transparent blue.

Additional Information

No.	SR-0166
Manufacturer (Brand)	Socorex
old neolab article no.	7-7436
EAN	4058072366858
Transport temperature	Room temperature
Color	blue
Material	Polycarbonate (PC)
Capacity	8
Type frame	Pipette rack
for pipettes / dispensers from	Socorex

